



***Evanthia:  
Leading Filler Technology***

Aesthetic Medicine made in Italy



We do not compromise when it comes to quality. Every stage of the production process is carefully monitored to ensure all parameters of the final product meet the conformity limits, with some being way below the accepted limits\*:

- Cross-linking agent residual:  
< 0,05 ppm
- Endotoxin content:  
≤ 0,5 Eu/ml
- Residual protein concentration:  
0,1%

\* Conformity parameters set by  
FDA and EU Regulation 2017/745.

Certifications of our manufacturing partner:  
UNI EN ISO 9001:2015 since 1995 with certificate No. 423.  
UNI CEI EN ISO 13485:2016 since 2008 with certificate No. 13790.  
ISO 22716:2007 since 2013 with certificate No. 18549.

Evanthia from Bioformula combines advanced technology with the safest manufacturing process – made in Italy, resulting in a safe and effective cross-linked hyaluronic acid dermal filler line for the face, that provides instant, long lasting, and natural looking results.

## *Evanthia 15*



Softens the appearance of fine lines.

Suggested treatments:

- Small periorcular wrinkles
- Small perioral wrinkles

## *Evanthia 20*



Redefines the face without losing its natural expressions.

Suggested treatments:

- Lip contours & enhancement
- Nasolabial folds
- Eyebrow lift
- Facial reshaping
- Marionette lines

## *Evanthia 25*



Enhances the face volume naturally and non-invasively.

Suggested treatments:

- Volumetric facial reshaping
- Lip enhancement
- Prejowl sulcus
- Deep nasolabial folds
- Temples
- Cheek augmentation
- Chin augmentation
- Rhinofiller
- Jawline contouring

***Evanthia:  
Suggested  
Treatments***



# *Evanthia:* *Technical* *Specifications\**

\* Please refer to package leaflet for full specifications.

\*\* Estimated average duration depends on various factors: patient's skin type, depth of wrinkles, injection area, and injected quantity. It is recommended to have the first two sessions within 4-6 months for the best result.

	<i><b>Evanthia 15 Basic</b></i>	<i><b>Evanthia 20 Medium</b></i>	<i><b>Evanthia 25 Superior</b></i>
Product code	40166	40167	40168
Composition	Reticulated hyaluronic acid of non-animal origin	Reticulated hyaluronic acid of non-animal origin	Reticulated hyaluronic acid of non-animal origin
Cross-linking agent	BDDE	BDDE	BDDE
HA concentration	15 mg/ml	20 mg/ml	25 mg/ml
Raw material molecular weight	2.400 – 4.000 kDa	2.400 – 4.000 kDa	2.400 – 4.000 kDa
Packaging	2 x 1ml	2 x 1ml	2 x 1ml
Blister	2 pre-filled, graduated, disposable, sterile glass syringes with Luer-Lok™ adapter	2 pre-filled, graduated, disposable, sterile glass syringes with Luer-Lok™ adapter	2 pre-filled, graduated, disposable, sterile glass syringes with Luer-Lok™ adapter
Needles	4 sterile disposable needles 30G 1/2" (CE 0123)	4 sterile disposable needles 27G 1/2" (CE 0123)	4 sterile disposable needles 27G 1/2" (CE 0123)
Injection depth	Superficial derma	Medium derma	Deep derma
Estimated duration of results**	Up to 6 months	Up to 8 months	Up to 10 months



[www.bioformula.it](http://www.bioformula.it)

Bioformula injectable solutions are Class III  
Medical Devices (CE 0373). For medical use only.  
Bioformula S.R.L. Viale delle Industrie 8  
20865 Usmate-Velate (MB) Italy +39 039 631 751

