


Plenhyage XL from Bioformula is a polynucleotidebased skin injectable that improves skin elasticity and is effective against aging, photoaging and for treating fibrous tissues. Its hydrating and anti-free radical activity help create the ideal environment for fibroblast growth, thereby restoring tissue elasticity and promoting cell turnover.

Plenhyage XL is the medical aesthetic treatment that naturally restores the optimal conditions of the skin for both face and body.



Plenhyage XL: Suggested Treatments



Plenhyage XL: Technical Specifications*

* Please refer to the package leaflet for full specifications.

	Plenhyage XL Medium	Plenhyage XL Strong
Product Code	40208	40209
Composition	Polymerised Polynucleotides	Polymerised Polynucleotides
PN Concentration	20 mg/ml	25 mg/ml
Origin of Raw Material	Wild Salmon	Wild Salmon
Molecular Weight of Raw Material	1000 kDa	1000 kDa
Packaging	1 x 2 ml	1 x 2 ml
Blister	1 pre-filled, graduated, disposable, sterile glass syringe with Luer-Lok™ adapter	1 pre-filled, graduated, disposable, sterile glass syringe with Luer-Lok™ adapter
Needles	2 sterile disposable needles 30G 1/2" (CE 0123)	2 sterile disposable needles 30G 1/2" (CE 0123)
Injection Depth	Derma	Derma



Origin of the Plenhylage XL Raw Material

Plenhyage XL is biocompatible and absorbable.

It features a premium, highly purified and safe raw material from Salmon sperm DNA that minimizes immunological side effects.





www.bioformula.it

Bioformula injectable solutions are Class III Medical Devices (CE 0373). For medical use only. Bioformula S.R.L. Viale delle Industrie 8 20865 Usmate-Velate (MB) Italy +39 039 631 751



